


Kii[®]
Abdominal Access Systems

Kii Fios®

First Entry Access System

THE ONLY INSUFFLATING OPTICAL PORT




Kii Fios
Access System

Insufflation
Needle

5mm
Access System


12mm
Access System

3mm



Minimal Fascial Defect

True nonbladed access allows for reapproximation of tissue layers with minimal fascial defect once the trocar is removed.




Minimal Penetration with Rapid Insufflation

Only 3mm of trocar penetration is needed to initiate insufflation.

Kii

Balloon Blunt-Tip Access System

INNOVATIVE BLUNT-TIP TECHNOLOGY



Now in the 5mm size – the smallest Hasson trocar available!

Suture Ties

Maximize instrument reach and facilitate port-site closure.

Non-Latex Balloon

Provides stable anchoring throughout lengthy procedures while maximizing working space and visibility.

Gel Cone

Enhances instrument articulation.


Optical Blunt-Tip Obturator

Facilitates direct visualization during trocar placement.

Kii


Advanced Fixation Technology

SUPERIOR RETENTION



Retention Disk

Reduces migration of the trocar sleeve into the abdominal cavity.



Non-Latex Balloon

Provides unmatched retention.
Maximizes working space and visibility.
Provides pressure on the incision site, potentially reducing port-site bleeding.
Maintains stability throughout procedures.

Kii

Low-Profile Access System

SLIM, LIGHTWEIGHT DESIGN



Low-Profile Seal Housing

Optimizes external working space, minimizing interference with adjacent ports.



Low-Profile Kii Seal



Standard 5mm Kii Seal

Z-Thread Sleeve with Extended Threads

Maintains retention at shorter depths, maximizing internal working space and visibility.

Kii Low-Profile Access System

Reorder No.		Size	
<i>Optical</i>	<i>Bladed</i>	<i>Sleeve</i>	
CTR21	CTB21	CTS28	5 x 55mm
CTR23	CTB23	N/A	5 x 55mm, dual pack
CTR22	CTB22	CTS29	5 x 100mm
CTR24	CTB24	N/A	5 x 100mm, dual pack

Kii

Shielded Bladed Access System

SMOOTH, CONTROLLED ENTRY

Intuitive Switch Design

Provides a clear indication of shield position.



Deactivated



Activated

Kii Shielded Bladed Access System


Reorder No.		Size
	<i>Advanced Fixation</i>	<i>Z-Thread</i>
N/A	CTB05	5 x 55mm
CFB03	CTB03	5 x 100mm
CFB01	CTB01	5 x 150mm
N/A	COQ61	8 x 100mm
CFB33	CTB33	11 x 100mm
CFB73	CTB73	12 x 100mm
CFB71	CTB71	12 x 150mm

Kii


Optical Access System

UNPARALLELED OPTICAL CLARITY


Tissue Reapproximation



Initial Withdrawal



Midway



Complete Withdrawal



Kii Optical Access System

Reorder No.		Size
	<i>Advanced Fixation</i>	<i>Z-Thread</i>
N/A	CTR05	5 x 55mm
CFR03	CTR03	5 x 100mm
CFR01	CTR01	5 x 150mm
N/A	COQ19	8 x 100mm
CFR33	CTR33	11 x 100mm
N/A	CTR34	11 x 100mm, handled
N/A	COB83	12 x 60mm
CFR73	CTR73	12 x 100mm
CFR74	CTR74	12 x 100mm, handled
CFR71	CTR71	12 x 150mm
N/A	CTR72	12 x 150mm, handled
N/A	COB37	15 x 100mm
CFR39	COB39	15 x 150mm



Kii Fios First Entry Access System

Reorder No.		Size
	<i>Advanced Fixation</i>	<i>Z-Thread</i>
CFF05	N/A	5 x 75mm
CFF03	CTF03	5 x 100mm
CFF04	CTF04	5 x 100mm, handled
CFF01	CTF01	5 x 150mm
CFF33	CTF33	11 x 100mm
CFF34	CTF34	11 x 100mm, handled
CFF73	CTF73	12 x 100mm
CFF74	CTF74	12 x 100mm, handled
CFF71	CTF71	12 x 150mm



Kii Balloon Blunt-Tip Access System

Reorder No.		Size
	<i>Blunt-Tip</i>	
COQ20		5 x 100mm
COR47		12 x 100mm
COR50		12 x 130mm



Kii Sleeve with Advanced Fixation

Reorder No.		Size
	<i>Advanced Fixation</i>	
CFS03		5 x 75mm
CFS02		5 x 100mm
CFS01		5 x 150mm
CFS12		11 x 100mm
CFS22		12 x 100mm
CFS21		12 x 150mm

Visit [appliedmedical.com/trocars](https://www.appliedmedical.com/trocars)
for more information.

© 2021 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.
254155-EN-USA-C